

Product datasheet for **TA310124**


PPAR gamma (PPARG) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant protein
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	peroxisome proliferator activated receptor gamma
Database Link:	NP_619725 Entrez Gene 25664 Rat Entrez Gene 5468 Human P37231
Synonyms:	CIMT1; GLM1; NR1C3; PPARG1; PPARG2; PPARgamma
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors
Protein Pathways:	Huntington's disease, Pathways in cancer, PPAR signaling pathway, Thyroid cancer



[View online »](#)

Product images:

Western blot analysis of rat testis lysates with anti-PPARG antibody (1:500).