

Product datasheet for **TA310133**


Thrombospondin 1 (THBS1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	WB: 1:50~200, ELISA: 1:40,000~80,000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Rabbit polyclonal THBS1 antibody was raised against a peptide from near the N terminal residues of human THBS1 protein.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Protein A purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	thrombospondin 1
Database Link:	NP_003237 Entrez Gene 21825 Mouse Entrez Gene 445442 Rat Entrez Gene 7057 Human P07996
Synonyms:	THBS; THBS-1; TSP; TSP-1; TSP1
Protein Families:	Druggable Genome
Protein Pathways:	Bladder cancer, ECM-receptor interaction, Focal adhesion, p53 signaling pathway, TGF-beta signaling pathway



[View online »](#)

Product images:

Immunohistochemical staining of formalin-fixed paraffin-embedded tonsil showing membrane staining with anti-THBS1 antibody diluted at 1:100.