

Product datasheet for **TA310139**


ABCB5 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Human ABCB5 (910-940 amino acid region)
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ATP binding cassette subfamily B member 5
Database Link:	NP_848654 Entrez Gene 340273 Human Q2M3G0
Synonyms:	ABCB5alpha; ABCB5beta; EST422562
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	ABC transporters



[View online »](#)

Product images:

Western blot analysis of Jurkat cells (Lane 1: 45ug lysate), Thp-1 (Lane 2: 60ug lysate) and ECL cells (Lane 3: 60ug) with anti-ABCB5 antibody diluted at 1:2000.