

Product datasheet for **TA310140**

ABCB9 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from human ABCB9.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ATP binding cassette subfamily B member 9
Database Link:	NP_062570 Entrez Gene 23457 Human Q9NP78
Synonyms:	EST122234; TAPL
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	ABC transporters, Lysosome



[View online »](#)

Product images:



Western blot analysis of THP1 cells (60ug lysate) with anti-ABCB9 antibody diluted at 1:2000. (Predicted and observed size: 80kDa)