

Product datasheet for **TA310141**


ABHD12 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Human ABHD12.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	abhydrolase domain containing 12
Database Link:	NP_001035937 Entrez Gene 76192 Mouse Entrez Gene 499913 Rat Entrez Gene 26090 Human Q8N2K0
Synonyms:	ABHD12A; BEM46L2; C20orf22; dj965G21.2; PHARC
Protein Families:	Protease, Transmembrane



[View online »](#)

Product images:

Western blot analysis of THP1 cells (Lanes 1 and 2: 60ug lysate) with anti-ABHD12 antibody diluted at 1:2000. (Predicted and observed size: 45kDa)