

Product datasheet for **TA310142**


ACOX2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Human ACOX2.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	acyl-CoA oxidase 2, branched chain
Database Link:	NP_003491 Entrez Gene 8309 Human Q99424
Synonyms:	BCOX; BRCACOX; BRCOX; THCCox
Protein Pathways:	Metabolic pathways, PPAR signaling pathway, Primary bile acid biosynthesis



[View online »](#)

Product images:

Western blot analysis of THP1 cells (Lane 1: 60ug lysate), ECL cells (Lane 2: 60ug lysate) and IM9 cells (Lane 3: 65ug lysate) with anti-ACOX2 antibody diluted at 1:2000. (Predicted and observed size: 77kDa)