

Product datasheet for **TA310146**

Acidic Calponin (CNN3) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:1000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Rabbit polyclonal CNN3 antibody was raised against a 16 amino acid peptide corresponding to the C terminus residues of human CNN3.
Formulation:	PBS, pH 7.4 with 0.02% sodium azide.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	calponin 3
Database Link:	NP_001830 Entrez Gene 1266 Human Q15417
Synonyms:	acidic; calponin; calponin 3; dj639P13.2.2 (acidic calponin 3)



[View online »](#)

Product images:



Western blot analysis of fetal heart lysates with anti-CNN3 antibody (1:200).