

Product datasheet for **TA311302**


USP6 Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1-3ug/ml
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Peptide with sequence KISPLHHLQMECSP, from the C Terminus of the protein sequence according to NP_004496.1.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20C. Minimize freezing and thawing.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ubiquitin specific peptidase 6
Database Link:	NP_004496 Entrez Gene 9098 Human P35125
Synonyms:	HRP1; Tre-2; TRE2; TRE17; TRESMCR; USP6-short
Protein Families:	Druggable Genome, Protease



[View online »](#)

Product images:

(2ug/ml) staining of Human Placenta lysate (35ug protein in RIPA buffer) with (B) and without (A) blocking with the immunizing peptide. Primary incubation was 1 hour. Detected by chemiluminescence.