

Product datasheet for **TA311321**


ABCD2 Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 0.5-2ug/ml
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Peptide with sequence C-PLSDTLAIK GKVID, from the internal region of the protein sequence according to NP_005155.1.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20C. Minimize freezing and thawing.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ATP binding cassette subfamily D member 2
Database Link:	NP_005155 Entrez Gene 225 Human Q9UBJ2
Synonyms:	ABC39; ALDL1; ALDR; ALDRP; hALDR
Protein Families:	Druggable Genome
Protein Pathways:	ABC transporters



[View online »](#)

Product images:

(1ug/ml) staining of Human Liver lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.