

Product datasheet for **TA311387**

IL17ra Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1-3ug/ml
Reactivity:	Rat (Expected from sequence similarity: Mouse)
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Peptide with sequence C-DQEKHGDDSKIN, from the internal region of the protein sequence according to NP_032385.1.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20C. Minimize freezing and thawing.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 17 receptor A
Database Link:	NP_001101353 Entrez Gene 16172 Mouse Entrez Gene 312679 Rat
Synonyms:	CD217; CDw217; hIL-17R; IL-17RA; IL17R; MGC10262



[View online »](#)

Product images:

TA311387 (1ug/ml) staining of Rat Kidney lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.