

## Product datasheet for **TA311442**

### TRPC4AP Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1-3ug/ml
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Peptide with sequence C-KLTQETYPNTYIFD , from the internal region of the protein sequence according to NP_056453.1; NP_955400.1.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20C. Minimize freezing and thawing.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	transient receptor potential cation channel subfamily C member 4 associated protein
Database Link:	<a href="#">NP_056453</a> <a href="#">Entrez Gene 26133 Human</a> <a href="#">Q8TEL6</a>
Synonyms:	C20orf188; PPP1R158; TRRP4AP; TRUSS
Protein Families:	Druggable Genome



[View online »](#)

**Product images:**

(1ug/ml) staining of U937 lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.