

Product datasheet for **TA311549**


ABHD14B Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 0.3-1ug/ml
Reactivity:	Human (Expected from sequence similarity: Mouse, Rat, Dog, Pig, Cow)
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Peptide with sequence HPCYLDKPEEWHT, from the internal region of the protein sequence according to NP_116139.1.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20C. Minimize freezing and thawing.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	abhydrolase domain containing 14B
Database Link:	NP_116139 Entrez Gene 76491 Mouse Entrez Gene 300983 Rat Entrez Gene 484744 Dog Entrez Gene 84836 Human Q96IU4
Synonyms:	CIB; HEL-S-299



[View online »](#)

Product images:

(0.3ug/ml) staining of Liver lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.