

Product datasheet for **TA311588**

EYA1 Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1-3ug/ml
Reactivity:	Human (Expected from sequence similarity: Mouse, Rat, Dog, Cow)
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Peptide with sequence C-TDPTAEYSTIHSP, from the internal region of the protein sequence according to NP_742057.1; NP_000494.2; NP_742056.1.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20C. Minimize freezing and thawing.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	EYA transcriptional coactivator and phosphatase 1
Database Link:	NP_742057 Entrez Gene 14048 Mouse Entrez Gene 502935 Rat Entrez Gene 477910 Dog Entrez Gene 2138 Human Q99502
Synonyms:	BOP; BOR; BOS1; OFC1
Protein Families:	Druggable Genome, Phosphatase, Transcription Factors



[View online »](#)

Product images:

(2ug/ml) staining of HEK293 lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.