

## Product datasheet for **TA311604**


### NPAS4 Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 0.3-1ug/ml
Reactivity:	Human, Mouse, Rat (Expected from sequence similarity: Dog, Cow)
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Peptide with sequence CRFNTSKSLRRQS, from the internal region of the protein sequence according to NP_849195.2.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20C. Minimize freezing and thawing.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	neuronal PAS domain protein 4
Database Link:	<a href="#">NP_849195</a> <a href="#">Entrez Gene 225872 Mouse</a> <a href="#">Entrez Gene 266734 Rat</a> <a href="#">Entrez Gene 483711 Dog</a> <a href="#">Entrez Gene 266743 Human</a> <a href="#">Q8IUM7</a>
Synonyms:	bHLHe79; Le-PAS; NXF; PASD10
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

**Product images:**

(0.3ug/ml) staining of Human Frontal Cortex lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.