

## Product datasheet for **TA311664**

### Mmp23 Rabbit Polyclonal Antibody


#### Product data:

|                       |  |
|-----------------------|--|
| Product Type:         | Primary Antibodies   |
| Applications:         | WB   |
| Recommended Dilution: | WB: 1:500~1:3000, ELISA: 1:10000   |
| Reactivity:           | Human, Mouse, Rat  |
| Host:                 | Rabbit   |
| Isotype:              | IgG  |
| Clonality:            | Polyclonal   |
| Immunogen:            | The antiserum was produced against synthesized peptide derived from human MMP-23 antibody.   |
| Formulation:          | Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Concentration:        | lot specific   |
| Purification:         | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                |
| Conjugation:          | Unconjugated   |
| Storage:              | Store at -20°C as received.  |
| Stability:            | Stable for 12 months from date of receipt.   |
| Gene Name:            | matrix metalloproteinase 23  |
| Database Link:        | <a href="#">NP_036115</a><br><a href="#">Entrez Gene 94339 Rat</a> <a href="#">Entrez Gene 26561 Mouse</a><br><a href="#">O88676</a> |
| Synonyms:             | CA-MMP   |
| Note:                 | MMP-23 antibody detects endogenous levels of total MMP-23 protein.   |



[View online »](#)

## Product images:



Western blot analysis of extracts from SK-OV3 cells, using MMP-23 antibody (#TA311664). The lane on the right is treated with the synthesized peptide.