

Product datasheet for **TA311704**

RPL26L1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from C-terminal of human RPL26L.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ribosomal protein L26 like 1
Database Link:	NP_057177 Entrez Gene 51121 Human Q9UNX3
Synonyms:	RPL26P1
Note:	RPL26L antibody detects endogenous levels of total RPL26L protein.
Protein Pathways:	Ribosome



[View online »](#)

Product images:



Western blot analysis of extracts from HUVEC cells and MCF-7 cells, using RPL26L antibody. The lane on the right is treated with the synthesized peptide.