

Product datasheet for **TA311762**

NCF1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:20000
Reactivity:	Human
Modifications:	Phospho-specific
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized phosphopeptide derived from human Neutrophil Cytosol Factor 1 around the phosphorylation site of serine 304 (R-S-SP-I-R).
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	neutrophil cytosolic factor 1
Database Link:	NP_000256 Entrez Gene 653361 Human P14598
Synonyms:	NCF1A; NOXO2; p47phox; SH3PXD1A
Note:	Neutrophil Cytosol Factor 1 (Phospho-Ser304) antibody detects endogenous levels of Neutrophil Cytosol Factor 1 only when phosphorylated at serine 304.



[View online »](#)

Protein Pathways:

Chemokine signaling pathway, Fc gamma R-mediated phagocytosis, Leukocyte transendothelial migration, Pathogenic Escherichia coli infection, Regulation of actin cytoskeleton