

Product datasheet for TA312003

MMP7 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 1:500~1:3000, IHC: 1:50~1:100, ELISA: 1:5000

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized peptide derived from human MMP-7.

Formulation: Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium

azide and 50% glycerol.

Concentration: lot specific

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: matrix metallopeptidase 7

Database Link: NP 002414

Entrez Gene 17393 MouseEntrez Gene 25335 RatEntrez Gene 4316 Human

P09237

Synonyms: MMP-7; MPSL1; PUMP-1

Note: MMP-7 antibody detects endogenous levels of total MMP-7 protein.

Protein Families: Druggable Genome, Protease

Protein Pathways: Wnt signaling pathway




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Product images:



Western blot analysis of extracts from COS7 cells, using MMP-7 antibody (#TA312003). The lane on the right is treated with the synthesized peptide.



Western blot analysis of extracts from COS7 cells (Lane 2) and HuvEc cells (Lane 3), using MMP-7 Antibody. The lane on the left is treated with systhesized peptide.



Immunohistochemical analysis of paraffinembedded human colon carcinoma tissue using MMP-7 antibody (#TA312003). The picture on the right is treated with the synthesized peptide.