

Product datasheet for **TA312390**

RAD18 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human RAD18.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	RAD18, E3 ubiquitin protein ligase
Database Link:	NP_064550 Entrez Gene 362412 Rat Entrez Gene 56852 Human Q9NS91
Synonyms:	RNF73
Note:	RAD18 antibody detects endogenous levels of total RAD18 protein.
Protein Families:	Druggable Genome



[View online »](#)

Product images:



Western blot analysis of extracts from LOVO cells, HepG2 cells and K562 cells, using RAD18 antibody. The lane on the right is treated with the synthesized peptide.