

Product datasheet for TA312587

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

STAT1 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 1:500~1:3000, IHC: 1:50~1:100, ELISA: 1:10000, IP: Various Dilution

Reactivity: Human, Mouse, Rat **Modifications:** Phospho-specific

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized phosphopeptide derived from human

STAT1 around the phosphorylation site of tyrosine 701 (T-G-YP-I-K).

Formulation: Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium

azide and 50% glycerol.

Concentration: lot specific

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: signal transducer and activator of transcription 1

Database Link: NP 009330

Entrez Gene 20846 MouseEntrez Gene 25124 RatEntrez Gene 6772 Human

P42224

Synonyms: CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91

Note: STAT1 (phospho-Tyr701) antibody detects endogenous levels of STAT1 only when

phosphorylated at tyrosine 701.

Protein Families: Druggable Genome, Transcription Factors







Protein Pathways:

Chemokine signaling pathway, Jak-STAT signaling pathway, Pancreatic cancer, Pathways in cancer, Toll-like receptor signaling pathway