

Product datasheet for **TA312663**


CHST6 Rabbit Polyclonal Antibody

Product data:

| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | IF |
| Recommended Dilution: | IF: 1:100~1:500, ELISA: 1:20000 |
| Reactivity: | Human |
| Host: | Rabbit |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Immunogen: | The antiserum was produced against synthesized peptide derived from internal of human CHST6. |
| Formulation: | Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Concentration: | lot specific |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Gene Name: | carbohydrate sulfotransferase 6 |
| Database Link: | NP_067628 Entrez Gene 4166 Human Q9GZX3 |
| Synonyms: | MCDC1 |
| Note: | CHST6 antibody detects endogenous levels of total CHST6 protein. |
| Protein Pathways: | Keratan sulfate biosynthesis, Metabolic pathways |



[View online »](#)

Product images:

Immunofluorescence analysis of A549 cells, using CHST6 antibody. The picture on the right is treated with the synthesized peptide.