

## Product datasheet for **TA312690**

### Septin 3 (SEPT3) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:20000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human SEPT3.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	septin 3
Database Link:	<a href="#">NP_061979</a> <a href="#">Entrez Gene 24050 Mouse</a> <a href="#">Entrez Gene 56003 Rat</a> <a href="#">Entrez Gene 55964 Human</a> <a href="#">Q9UH03</a>
Synonyms:	bK250D10.3; SEP3
Note:	SEPT3 antibody detects endogenous levels of total SEPT3 protein.



[View online »](#)

Product images:



Western blot analysis of extracts from 293 cells, using SEPT3 antibody. The lane on the right is treated with the synthesized peptide.