

Product datasheet for **TA312717**

PASK Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human ATPG.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	PAS domain containing serine/threonine kinase
Database Link:	NP_055963 Entrez Gene 301617 Rat Entrez Gene 23178 Human Q96RG2
Synonyms:	PASKIN; STK37
Note:	ATPG antibody detects endogenous levels of total ATPG protein.
Protein Families:	Druggable Genome, Protein Kinase, Stem cell - Pluripotency



[View online »](#)

Product images:



Western blot analysis of extracts from HeLa cells and Jurkat cells, using ATPG antibody. The lane on the right is treated with the synthesized peptide.