

## Product datasheet for **TA312772**

### **Kallikrein 8 (KLK8) Rabbit Polyclonal Antibody**


#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	WB
<b>Recommended Dilution:</b>	WB: 1:500~1:3000, ELISA: 1:10000
<b>Reactivity:</b>	Human
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Neuropsin.
<b>Formulation:</b>	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Concentration:</b>	lot specific
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	kallikrein related peptidase 8
<b>Database Link:</b>	<a href="#">NP_009127</a> <a href="#">Entrez Gene 11202 Human</a> <a href="#">O60259</a>
<b>Synonyms:</b>	HNP; NP; NRPN; PRSS19; TADG14
<b>Note:</b>	Neuropsin (Cleaved-Val33) antibody detects endogenous levels of fragment of activated Neuropsin resulting from cleavage adjacent to Val33.
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane



[View online »](#)

## Product images:



Western blot analysis of extracts from Jurkat cells, treated with etoposide (25uM, 24hours), using Neuropsin (Cleaved-Val33) antibody. The lane on the right is treated with the synthesized peptide.