

Product datasheet for TA312815

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PKC theta (PRKCQ) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB: 1:500~1:3000, IHC: 1:50~1:100, ELISA: 1:20000

Reactivity: Human, Mouse, Rat **Modifications:** Phospho-specific

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized phosphopeptide derived from human PKC?

around the phosphorylation site of serine 676 (R-L-SP-F-A).

Formulation: Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium

azide and 50% glycerol.

Concentration: lot specific

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: protein kinase C theta

Database Link: NP 006248

Entrez Gene 18761 MouseEntrez Gene 85420 RatEntrez Gene 5588 Human

Q04759

Synonyms: nPKC-theta; PRKCT

Note: PKC. (phospho-Ser676) antibody detects endogenous levels of PKC. only when

phosphorylated at serine 676.

Protein Families: Druggable Genome, Protein Kinase, Transcription Factors







Protein Pathways:

Adipocytokine signaling pathway, T cell receptor signaling pathway, Tight junction, Vascular smooth muscle contraction