

## Product datasheet for **TA312901**

### JUND Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1:500~1:3000, IHC: 1:50~1:100, ELISA: 1:20000, IP: Various Dilution
Reactivity:	Human, Mouse, Rat
Modifications:	Phospho-specific
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized phosphopeptide derived from human JunD around the phosphorylation site of serine 255 (G-E-SP-P-P).
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	JunD proto-oncogene, AP-1 transcription factor subunit
Database Link:	<a href="#">NP_005345</a> <a href="#">Entrez Gene 16478 Mouse</a> <a href="#">Entrez Gene 24518 Rat</a> <a href="#">Entrez Gene 3727 Human</a> <a href="#">P17535</a>
Synonyms:	AP-1
Note:	JunD (phospho-Ser255) antibody detects endogenous levels of JunD only when phosphorylated at serine 255.
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	MAPK signaling pathway



[View online »](#)