

Product datasheet for TA313026

TAK1 (MAP3K7) Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1:500~1:3000, IHC: 1:50~1:100, ELISA: 1:1000
Reactivity:	Human, Mouse, Rat
Modifications:	Phospho-specific
Host:	Rabbit
lsotype:	lgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized phosphopeptide derived from human Stathmin around the phosphorylation site of serine 439 (R-R-SP-I-Q).
Formulation:	Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mitogen-activated protein kinase kinase kinase 7
Database Link:	<u>NP_003179</u> <u>Entrez Gene 26409 MouseEntrez Gene 313121 RatEntrez Gene 6885 Human</u> <u>O43318</u>
Synonyms:	MEKK7; TAK1; TGF1a
Note:	MAP3K7 (Phospho-Ser439) antibody detects endogenous levels of MAP3K7 only when phosphorylated at serine 439.
Protein Families:	Druggable Genome, Protein Kinase



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US Protein Pathways:Adherens junction, MAPK signaling pathway, NOD-like receptor signaling pathway, RIG-I-like
receptor signaling pathway, T cell receptor signaling pathway, Toll-like receptor signaling
pathway, Wnt signaling pathway

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US