

Product datasheet for TA313711

Cathepsin L (CTSL) Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

| Due du et Ture et | |
|-----------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Recommended Dilution: | WB: 1:500~1:3000, ELISA: 1:10000 |
| Reactivity: | Human |
| Host: | Rabbit |
| lsotype: | lgG |
| Clonality: | Polyclonal |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human CATL1. |
| Formulation: | Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Concentration: | lot specific |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Gene Name: | cathepsin L |
| Database Link: | <u>NP_666023</u> <u>Entrez Gene 1514 Human</u> <u>P07711</u> |
| Synonyms: | CATL; CTSL1; MEP |
| Note: | CATL1 (heavy chain, Cleaved-Thr288) antibody detects endogenous levels of fragment of activated CATL1 resulting from cleavage adjacent to Thr288. |
| Protein Families: | Druggable Genome, Protease |
| Protein Pathways: | Antigen processing and presentation, Lysosome |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:



Western blot analysis of extracts from HeLa cells, treated with etoposide (25uM, 1hour), using CATL1 (heavy chain, Cleaved-Thr288) antibody.The lane on the right is treated with the synthesized peptide.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US