

Product datasheet for **TA313804**

DAP5 (EIF4G2) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:20000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human EIF4G2.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	eukaryotic translation initiation factor 4 gamma 2
Database Link:	NP_001036024 Entrez Gene 13690 Mouse Entrez Gene 1982 Human P78344
Synonyms:	AAG1; DAP5; NAT1; P97
Note:	EIF4G2 antibody detects endogenous levels of total EIF4G2 protein.
Protein Families:	Transcription Factors
Protein Pathways:	Viral myocarditis



[View online »](#)

Product images:



Western blot analysis of extracts from A549 cells, using EIF4G2 antibody. The lane on the right is treated with the synthesized peptide.