

Product datasheet for **TA313870**

FER Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1:500~1:3000, IHC: 1:50~1:100, ELISA: 1:20000
Reactivity:	Human, Mouse
Modifications:	Phospho-specific
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized phosphopeptide derived from human FER around the phosphorylation site of tyrosine 402 (V-N-YP-E-E).
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	FER tyrosine kinase
Database Link:	NP_005237 Entrez Gene 14158 Mouse Entrez Gene 103921 Mouse Entrez Gene 2241 Human P16591
Synonyms:	FerT; Pe1Fe3; Pe1Fe6; Pe1Fe10; Pe1Fe13; PPP1R74; TYK3
Note:	FER (Phospho-Tyr402) antibody detects endogenous levels of FER only when phosphorylated at tyrosine 402.
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Adherens junction



[View online »](#)