

Product datasheet for **TA314069**

IL2 Receptor beta (IL2RB) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:10000
Reactivity:	Human, Mouse
Modifications:	Phospho-specific
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized phosphopeptide derived from human IL-2R β around the phosphorylation site of tyrosine 364 (Q-G-YP-F-F).
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 2 receptor subunit beta
Database Link:	NP_000869 Entrez Gene 16185 Mouse Entrez Gene 3560 Human P14784
Synonyms:	CD122; IL15RB; P70-75
Note:	IL-2R. (Phospho-Tyr364) antibody detects endogenous levels of IL-2R. only when phosphorylated at tyrosine 364.
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, Endocytosis, Jak-STAT signaling pathway



[View online »](#)