

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA314095

CD51 (ITGAV) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:10000
Reactivity:	Human
Host:	Rabbit
lsotype:	lgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from human ITGAV.
Formulation:	Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	integrin subunit alpha V
Database Link:	<u>NP_001138471</u> <u>Entrez Gene 3685 Human</u> <u>P06756</u>
Synonyms:	CD51; MSK8; VNRA; VTNR
Note:	ITGAV (heavy chain, Cleaved-Lys889) antibody detects endogenous levels of fragment of activated ITGAV resulting from cleavage adjacent to Lys889.
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

CD51 (ITGAV) Rabbit Polyclonal Antibody – TA314095

Protein Pathways:Arrhythmogenic right ventricular cardiomyopathy (ARVC), Cell adhesion molecules (CAMs),
Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hypertrophic
cardiomyopathy (HCM), Pathways in cancer, Regulation of actin cytoskeleton, Small cell lung
cancer

Product images:

A549 A5	549	
	117	
ITGAV (heavy chain,Cleaved-Lys889)	85	
	48	Western blot analysis of extracts from A549 cells, treated with etoposide (25uM, 1hour), using ITGAV (heavy chain, Cleaved-Lys889) antibody.The lane on the right is treated with the
	34	synthesized peptide.
	26	
	19	
	(kD)	

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US