

Product datasheet for **TA314205**

MMP10 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from human MMP10.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	matrix metalloproteinase 10
Database Link:	NP_002416 Entrez Gene 4319 Human P09238
Synonyms:	SL-2; STMY2
Note:	MMP10 (Cleaved-Phe99) antibody detects endogenous levels of fragment of activated MMP10 resulting from cleavage adjacent to Phe99.
Protein Families:	Druggable Genome, Protease



[View online »](#)

Product images:



Western blot analysis of extracts from HepG2 cells, using MMP10 (Cleaved-Phe99) antibody. The lane on the right is treated with the synthesized peptide.