

# Product datasheet for TA314354

## PML Protein (PML) Rabbit Polyclonal Antibody

## **Product data:**

#### **Product Type: Primary Antibodies Applications:** IF, IHC, WB **Recommended Dilution:** WB: 1:500~1:3000, IHC: 1:50~1:100, IF: 1:100~1:500, ELISA: 1:10000 Reactivity: Human Rabbit Host: Isotype: lgG **Clonality:** Polyclonal Immunogen: The antiserum was produced against synthesized peptide derived from human PML. Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium Formulation: azide and 50% glycerol. **Concentration:** lot specific **Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. Unconjugated **Conjugation:** Store at -20°C as received. Storage: Stable for 12 months from date of receipt. Stability: Gene Name: promyelocytic leukemia Database Link: NP 002666 Entrez Gene 5371 Human P29590 Synonyms: MYL; PP8675; RNF71; TRIM19 Note: PML antibody detects endogenous levels of total PML protein. **Protein Families:** Druggable Genome, Transcription Factors **Protein Pathways:** Acute myeloid leukemia, Pathways in cancer, Ubiquitin mediated proteolysis




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn


## **Product images:**



Western blot analysis of extracts from A549 cells, using PML antibody (#TA314354). The lane on the right is treated with the synthesized peptide.




Western blot analysis of extracts from JK cells (Lane 2), using PML Antibody. The lane on the left is treated with synthesized peptide.



Immunohistochemical analysis of paraffinembedded human tonsil tissue using PML antibody (#TA314354).The picture on the right is treated with the synthesized peptide.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US





Immunofluorescence analysis of A549 cells, using PML antibody (#TA314354).The picture on the right is treated with the synthesized peptide.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US