

Product datasheet for TA314498

RPS25 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

IHC **Applications:**

Recommended Dilution: IHC: 1:50~1:100, ELISA: 1:10000

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: lgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

RPS25.

Formulation: Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium

azide and 50% glycerol.

Concentration: lot specific

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Conjugation: Unconjugated

Store at -20°C as received. Storage:

Stability: Stable for 12 months from date of receipt.

Gene Name: ribosomal protein S25

Database Link: NP 001019

Entrez Gene 75617 MouseEntrez Gene 122799 RatEntrez Gene 6230 Human

P62851

Synonyms: S25

RPS25 antibody detects endogenous levels of total RPS25 protein. Note:

Protein Pathways: Ribosome




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Product images:



Immunohistochemistry analysis of paraffinembedded human breast carcinoma tissue using RPS25 antibody. The picture on the right is treated with the synthesized peptide.