

Product datasheet for TA314556

SOX22 (SOX12) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC

Recommended Dilution: IHC: 1:50~1:100, ELISA: 1:20000

Reactivity: Human, Mouse

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized peptide derived from internal of human

SOX12

Formulation: Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium

azide and 50% glycerol.

Concentration: lot specific

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: SRY-box 12

Database Link: NP 008874

Entrez Gene 20667 MouseEntrez Gene 6666 Human

O15370

Synonyms: SOX22

Note: SOX12 antibody detects endogenous levels of total SOX12 protein.

Protein Families: ES Cell Differentiation/IPS, Transcription Factors




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Product images:



Immunohistochemistry analysis of paraffinembedded human brain tissue, using SOX12 antibody. The picture on the right is treated with the synthesized peptide.