

Product datasheet for **TA314632**

TNFAIP1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human TNAP1.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	TNF alpha induced protein 1
Database Link:	NP_066960 Entrez Gene 21927 Mouse Entrez Gene 287543 Rat Entrez Gene 7126 Human Q13829
Synonyms:	B12; B61; BTBD34; EDP1; hBACURD2
Note:	TNAP1 antibody detects endogenous levels of total TNAP1 protein.
Protein Families:	Druggable Genome, Ion Channels: Other



[View online »](#)

Product images:



Western blot analysis of extracts from HUVEC cells treated with PMA (125ng/ml, 30mins), and HUVEC cells treated with serum (20%, 30mins), using TNAP1 antibody. The lane on the right is treated with the synthesized peptide.