

Product datasheet for **TA314799**


SERPINB7 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:20000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human SERPINB7.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	serpin family B member 7
Database Link:	NP_001035237 Entrez Gene 8710 Human Q75635
Synonyms:	MEGSIN; PPKN; TP55
Note:	SERPINB7 antibody detects endogenous levels of total SERPINB7 protein.
Protein Families:	Druggable Genome

[View online »](#)

Product images:



Western blot analysis of extracts from Jurkat cells, using SERPINB7 antibody. The lane on the right is treated with the synthesized peptide.