

## Product datasheet for **TA314808**

### **RIOK3 Rabbit Polyclonal Antibody**


#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	WB
<b>Recommended Dilution:</b>	WB: 1:500~1:3000, ELISA: 1:40000
<b>Reactivity:</b>	Human, Mouse
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from internal of human RIOK3.
<b>Formulation:</b>	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Concentration:</b>	lot specific
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	RIO kinase 3
<b>Database Link:</b>	<a href="#">NP_003822</a> <a href="#">Entrez Gene 66878 Mouse</a> <a href="#">Entrez Gene 8780 Human</a> <a href="#">O14730</a>
<b>Synonyms:</b>	SUDD
<b>Note:</b>	RIOK3 antibody detects endogenous levels of total RIOK3 protein.
<b>Protein Families:</b>	Druggable Genome, Protein Kinase



[View online »](#)

## Product images:



Western blot analysis of extracts from HepG2 cells, using RIOK3 antibody. The lane on the right is treated with the synthesized peptide.