

## Product datasheet for **TA314811**

### IL18R1 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:10000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human IL18R.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 18 receptor 1
Database Link:	<a href="#">NP_003846</a> <a href="#">Entrez Gene 8809 Human</a> <a href="#">Q13478</a>
Synonyms:	CD218a; CDw218a; IL-1Rrp; IL1RRP; IL18RA
Note:	IL18R antibody detects endogenous levels of total IL18R protein.
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction



[View online »](#)

## Product images:



Western blot analysis of extracts from Jurkat cells, using IL18R antibody. The lane on the right is treated with the synthesized peptide.