

Product datasheet for TA314836

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

WAVE 1 (WASF1) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IF, IHC, WB

Recommended Dilution: WB: 1:500~1:3000, IHC: 1:50~1:100, IF: 1:100~1:500, ELISA: 1:40000

Reactivity: Human, Mouse

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized non-phosphopeptide derived from human

WAVE1 around the phosphorylation site of tyrosine 125 (E-T-YP-D-V).

Formulation: Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium

azide and 50% glycerol.

Concentration: lot specific

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

epitope-specific immunogen.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: WAS protein family member 1

Database Link: NP 001020105

Entrez Gene 83767 MouseEntrez Gene 8936 Human

Q92558

Synonyms: SCAR1; WAVE; WAVE1

Note: WAVE1 (Ab-125) antibody detects endogenous levels of total WAVE1 protein.


Protein Families: Druggable Genome

Protein Pathways: Adherens junction, Fc gamma R-mediated phagocytosis, Regulation of actin cytoskeleton







Product images:



Western blot analysis of extracts from 293 cells, treated with insulin (0.01U/ml, 15mins), using WAVE1 (Ab-125) antibody. The lane on the right is treated with the synthesized peptide.



Immunohistochemistry analysis of paraffinembedded human brain tissue using WAVE1 (Ab-125) antibody. The picture on the right is treated with the synthesized peptide.



Immunofluorescence analysis of COS-7 cells, using NFAT3 (Ab-676) antibody. The picture on the right is treated with the synthesized peptide.