

## Product datasheet for **TA314862**


### **RABEP1 Rabbit Polyclonal Antibody**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	IHC
<b>Recommended Dilution:</b>	IHC: 1:50~1:100, ELISA: 1:40000
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from internal of human RABEP1.
<b>Formulation:</b>	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Concentration:</b>	lot specific
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	rabaptin, RAB GTPase binding effector protein 1
<b>Database Link:</b>	<a href="#">NP_004694</a> <a href="#">Entrez Gene 54189 Mouse</a> <a href="#">Entrez Gene 54190 Rat</a> <a href="#">Entrez Gene 9135 Human</a> <a href="#">Q15276</a>
<b>Synonyms:</b>	RAB5EP; RABPT5
<b>Note:</b>	RABEP1 antibody detects endogenous levels of total RABEP1 protein.
<b>Protein Pathways:</b>	Endocytosis



[View online »](#)

**Product images:**

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using RABEP1 antibody. The picture on the right is treated with the synthesized peptide.