

Product datasheet for **TA315058**

Follistatin (FST) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:5000
Reactivity:	Mouse, Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human FST.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	follistatin
Database Link:	NP_006341 Entrez Gene 14313 Mouse Entrez Gene 24373 Rat Entrez Gene 10468 Human P19883
Synonyms:	FS
Note:	FST antibody detects endogenous levels of total FST protein.
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	TGF-beta signaling pathway



[View online »](#)

Product images:



Western blot analysis of extracts from HepG2, K562, COLO205, Jurkat and HeLa cells, using FST antibody. The lane on the right is treated with the synthesized peptide.