

Product datasheet for **TA315077**

TCP1 eta (CCT7) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:5000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human CCT7.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	chaperonin containing TCP1 subunit 7
Database Link:	NP_001009570 Entrez Gene 12468 MouseEntrez Gene 297406 RatEntrez Gene 10574 Human Q99832
Synonyms:	CCTETA; CCTH; NIP7-1; TCP1ETA
Note:	CCT7 antibody detects endogenous levels of total CCT7 protein.



[View online »](#)

Product images:



Western blot analysis of extracts from HeLa cells, using CCT7 antibody. The lane on the right is treated with the synthesized peptide.