

## Product datasheet for **TA315110**

### UTP14A Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:20000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human UTP14A.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	UTP14A small subunit processome component
Database Link:	<a href="#">NP_001159693</a> <a href="#">Entrez Gene 72554 Mouse</a> <a href="#">Entrez Gene 10813 Human</a> <a href="#">Q9BVJ6</a>
Synonyms:	dJ537K23.3; NYCO16; SDCCAG16
Note:	UTP14A antibody detects endogenous levels of total UTP14A protein.



[View online »](#)

## Product images:



Western blot analysis of extracts from HeLa cells, using UTP14A antibody. The lane on the right is treated with the synthesized peptide.