

Product datasheet for **TA315291**

ABHD14A Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Mouse, Rat, Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human ABHD14A.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	abhydrolase domain containing 14A
Database Link:	NP_056222 Entrez Gene 68644 MouseEntrez Gene 300982 RatEntrez Gene 25864 Human Q9BUJ0
Synonyms:	DORZ1
Note:	ABHD14A antibody detects endogenous levels of total ABHD14A protein.



[View online »](#)

Product images:



Western blot analysis of extracts from RAW264.7 cells, using ABHD14A antibody. The lane on the right is treated with the synthesized peptide.