

## Product datasheet for **TA315307**


### DECR2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:10000
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human DECR2.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	2,4-dienoyl-CoA reductase 2, peroxisomal
Database Link:	<a href="#">NP_065715</a> <a href="#">Entrez Gene 26378 MouseEntrez Gene 64461 RatEntrez Gene 26063 Human Q9NUI1</a>
Synonyms:	PDCR; SDR17C1
Note:	DECR2 antibody detects endogenous levels of total DECR2 protein.
Protein Families:	Druggable Genome



[View online »](#)

**Product images:**

Western blot analysis of extracts from LOVO cells, using DECR2 antibody. The lane on the right is treated with the synthesized peptide.