

Product datasheet for **TA315357**

MLH3 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:40000
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human MLH3.
Formulation:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mutL homolog 3
Database Link:	NP_001035197 Entrez Gene 27030 Human Q9UHC1
Synonyms:	HNPCC7
Note:	MLH3 antibody detects endogenous levels of total MLH3 protein.
Protein Families:	Druggable Genome
Protein Pathways:	Mismatch repair



[View online »](#)

Product images:



Western blot analysis of extracts from K562 cells and HT-29 cells, using MLH3 antibody. The lane on the right is treated with the synthesized peptide.