

## Product datasheet for **TA315371**

### GPR82 Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:500~1:3000, ELISA: 1:10000
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human GPR82.
Formulation:	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	G protein-coupled receptor 82
Database Link:	<a href="#">NP_543007</a> <a href="#">Entrez Gene 319200 Mouse</a> <a href="#">Entrez Gene 27197 Human</a> <a href="#">Q96P67</a>
Synonyms:	G protein-coupled receptor 82; G protein-coupled receptor GPR82; OTTMUSP00000018382
Note:	GPR82 antibody detects endogenous levels of total GPR82 protein.
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)

## Product images:



Western blot analysis of extracts from K562 cells, A549 cells and HT-29 cells, using GPR82 antibody. The lane on the right is treated with the synthesized peptide.