

Product datasheet for TA315462

CYP39A1 Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn


Product Type:	Primary Antibodies
Applications:	IF, WB
Recommended Dilution:	WB: 1:500~1:3000, IF: 1:100~1:500, ELISA: 1:10000
Reactivity:	Human, Mouse
Host:	Rabbit
lsotype:	IgG
Clonality:	Polyclonal
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human Cytochrome P450 39A1.
Formulation:	Phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration:	lot specific
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	cytochrome P450 family 39 subfamily A member 1
Database Link:	<u>NP_057677</u> <u>Entrez Gene 56050 MouseEntrez Gene 51302 Human</u> <u>Q9NYL5</u>
Synonyms:	39a1 (oxysterol 7alpha-hydroxylase); cytochrome P450; family 39; oxysterol 7-alpha- hydroxylase; oxysterol 7alpha-hydroxylase; polypeptide 1; subfamily a
Note:	Cytochrome P450 39A1 antibody detects endogenous levels of total Cytochrome P450 39A1 protein.
Protein Families:	Druggable Genome, P450, Transmembrane
Protein Pathways:	Primary bile acid biosynthesis




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Product images:



Western blot analysis of extracts from LOVO cells and HUVEC cells, using Cytochrome P450 39A1 antibody.The lane on the right is treated with the synthesized peptide.



Immunofluorescence analysis of A549 cells, using Cytochrome P450 39A1 antibody.The picture on the right is treated with the synthesized peptide.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US